

3 June, 2022

Re: ÄKTA™ ready Flow Kit Storage Conditions

To Whom It May Concern,

The recommended storage temperature for standard ÄKTA ready flow kits is provided in Section 8.3 of the Operating Instructions 28960345 and specified as $> +5^{\circ}\text{C}$. This recommendation also applies to all modified ÄKTA ready flow kits as well, including the two listed in the below table. Extended storage below the recommended $+5^{\circ}\text{C}$ could lead to brittleness or cracking of the plastic connectors. However, the operating temperature of ÄKTA ready flow kits is $+2^{\circ}\text{C}$ to $+40^{\circ}\text{C}$. If the kits are allowed to acclimate to a warmer temperature before being used this would reduce the risk of damage to the kit during setup and handling.

Description	Part Number	Operating Temperature
High Flow Kit F, Modified, ÄKTA ready	29477427	$+2^{\circ}\text{C}$ to $+40^{\circ}\text{C}$
High Flow Gradient C, Modified, ÄKTA ready	29184612	$+2^{\circ}\text{C}$ to $+40^{\circ}\text{C}$

Operating Instructions access link:

<https://cdn.cytivalifesciences.com/api/public/content/digi-14501-pdf>

Contact Information

This letter is provided for informational purposes only and does not affect or create any rights or obligations in regard to either you or Cytiva.

If you have any questions, please do not hesitate to contact your Cytiva representative.

Best regards,



Mike Toner
Product Manager
michaeltoner@cytiva.com

Certificate of Quality

Product: ÄKTA™ ready Gradient Flow Section With Inlets

Lot Number: 17242818

Product Article Number: 29-1846-12

Date of Manufacture: 20210126

Product Description: High Flow Gradient C, Modified, ÄKTA™ ready

Expiration Date: 20230126

Product Release Criteria

We hereby certify that the defined product has been manufactured to meet its specifications and have been verified to meet predetermined Critical to Quality attributes and the flow kit has been formally released for delivery.

Description	Limits
Autoclave Pump tubing	121°C > 15 min
Gamma Irradiation Inlets	25.0 – 40.0 kGy

Regulatory Conformance

Bio-safety

All polymeric materials in contact with the process fluid comply with United States Pharmacopeia (USP) <88> Biological Reactivity Test Class VI and are free from animal derived ingredients or in compliance with EMA/410/01.

Country of Origin: USA

This product is manufactured in compliance with our ISO 9001 certified quality management system.

Issued by Cytiva Westborough Quality Assurance

This document has been electronically produced and is valid without a signature.

ÄKTA is a trademark of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate.



Global Life Sciences Solutions USA LLC
14 Walkup Drive
Westborough, MA 01581 USA

Certificate of Quality

This product is manufactured in compliance with our ISO 9001 certified quality management system.

Issued by Cytiva Westborough Quality Assurance

This document has been electronically produced and is valid without a signature.

Product: AKTA ready Low Flow Kit
Lot Number: 18356721
Product Article Number: 28 9301 82
Date of Manufacture: 20240315
Product Description: Low Flow Kit, AKTA ready
Expiration Date: 20260315

Product Release Criteria

We hereby certify that the defined product has been manufactured to meet its specifications and have been verified to meet predetermined Critical to Quality attributes.

Description	Test Method	Result
Autoclave - Pump tubing	121C > 15 min	Conforms
Gamma Irradiation - Inlets	25.0 - 40.0 kGy	Conforms
Flow Rate Test	Per SOP-QC-042	Pass
Pressure Integrity	Max 0.5 bar, 5 min hold	Pass
Visual Inspection	No defects visible	Pass
Package Integrity	Sealed, no damage	Pass

Regulatory Conformance

Bio-safety: All polymeric materials in contact with the process fluid comply with United States Pharmacopeia (USP) <88> Biological Reactivity Test Class VI.

AKTA is a trademark of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate.